



## Clinton 8890 with Step Stool - Fully Refurbished

This versatile table has storage access from either side. Plus, its lower treatment table height could eliminate the need for a step stool, and makes it ideal for family practice, clinic or wherever a cost effective treatment/exam table is needed. It is available in your choice of 20 colors of upholstery with color coordinated laminate doors and drawers on an easy-clean laminate base. Includes a built-in, pull out step stool for easy patient access. Step stool slides back in when not in use.

### Features:

- Premium, vinyl upholstery
- Handy, pull-out, welded steel, step stool
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Optional pull-out basin
- Easy care all laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides
- 2 drawers with all steel slides
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included



### Specs:

- Dimensions: 58"/72"(L) x 28.5"(W) x 31"(H)
- Step Stool Dimensions: 16.875"(L) x 19.5"(W) x 7.75"(H)
- 400 lbs. (181.43 kg) load capacity under normal use
- 2" firm, foam padding (5 cm)

### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.