



Clinton Imaging Table- 3-Section Top/Drop Window - Fully Refurbished

The Clinton Imaging Table with Three-Section Top and Drop Window is a durable and reliable table that meets clinical standards and ensures patient comfort and satisfaction. The Imaging Table with Three-Section Top and Drop Window features a three-piece top with a pneumatic adjustable backrest. A manual footrest and a window drop are incorporated into each chair. The footrests are adjustable from either side and raise up to 30 degrees from a flat position. This table is constructed with an unobstructed base design.

Features:

- Three-piece top with pneumatic adjustable backrest, manual footrest and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails & casters
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- 34" wide model available
- 550 lbs load capacity under normal use

Specs:

- Dimensions: 76"(L) x 27"(W) x 21.5" – 37.5"(H)
- Backrest: 38"(L) x 27"(W)
- Footrest: 12"(L) x 27"(W)
- Window: 10"(L) x 10"(W)



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.