



Clinton Industries 3 Section Power Ultrasound- Fully Refurbished

Clinton Power Ultrasound Trendelenburg Tables are capable of trendelenburg and reverse trendelenburg positions within plus or minus of 15° for a wide range of uses. The Power Exam Tables feature a three-piece top which includes a pneumatic adjustable backrest and legrest. The backrest on the exam table can be adjusted from 30° to 60° on either side of the table. Easy-to-use hand controls operate the dual lift columns allowing the Ultrasound Table to be adjusted from 21.5 inches to 37.5 inches.

Features:

- Three-piece top with pneumatic adjustable backrest and manual footrest
- Footrest is adjustable from either side and raises up to 30° from flat position
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- Optional side rails & casters available
- 34" wide model available
- 550 lbs load capacity under normal use



Specs:

Dimensions: 76"(L) x 27"(W) x 21.5" – 37.5"(H)

Backrest: 38"(L) x 27"(W)

Footrest: 13"(L) x 27"(W)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.