



Clinton Industries Shrouded 34" Bariatric Power Table - Fully Refurbished

Clinton Industries Shrouded Extra Wide Bariatric Power Table With Adjustable Backrest And Drop Section, is specially designed with dual-cylinder pneumatic backrest, a drop section with a nose hole, an easy-clean ABS plastic enclosed lift, and an electric power height adjuster with a foot control pedal. The ultra-strong uni-frame base and top construction has a heavy duty all welded lift mechanism, and a durable gray powder-coated steel frame.

Features:

- Extra wide, bariatric model with choice of two widths
- Back rest is adjustable from either side and raises 80 degrees from flat
- Drop Section is adjustable from either side and raises 80 degrees from flat
- Standard drop section has nose hole, (can be ordered without a nose hole)
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- Low height may help achieve ADA tax credit
- Customize with options



Specs:

Dimension: 76"(L) x 34"(W) x 21" – 34"(H)

Backrest: 38"(L) x 34"(W)

Drop Section: 13"(L) x 34"(W)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.