



Clinton Industries Shrouded Power Table - Fully Refurbished

Clinton provides plenty of smooth dependable power to our line of Shrouded Power Tables. Engineered and built tough in the USA by Clinton Industries, the tables provide durability, convenience and versatility for a wide variety of applications. And not only do Clinton Shrouded Power Tables perform at superior levels, but they look great too! This table comes with a pneumatic backrest that is adjustable from either side and raises 80° from flat.

Features:

- Easy-clean, ABS plastic enclosed lift
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- 450 lbs. (204.11 kg) load capacity under normal use
- Dimensions: 72"(L) x 27"(W) x 21" – 34"(H)
- Backrest: 38"(L) x 27"(W)



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.