



Clinton Industries Shrouded Power Table - Fully Refurbished

The Medical Examination Table features an adjustable, pneumatic backrest that can be adjusted up to 80 degrees from a flat position. The Exam Table comes with retractable, 2-position, length adjustable stirrups. The 81360 Power Table has an enclosed scissor lift mechanism. The Clinton Electric Table offers a drop section that is also adjustable from both sides; and raises up to 30 degrees from the flat position.

Features:

- Pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- Footrest is adjustable from either side and raises 80 degrees from flat
- Easy-clean, ABS plastic enclosed lift
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top



Specs:

- Dimensions: 76"(L) x 27"(W) x 21" – 34"(H)
- Backrest: 38"(L) x 27"(W)
- Footrest: 13"(L) x 19"(W)
- 2" (5 cm) ultra-firm padding
- 450 lbs. (204.11 kg) load capacity under normal use

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.