



Imaging Table with Three-Section Top and Drop Window

Clinton offers a broad range of medical tables designed to meet the needs and preferences of physicians, hospitals and facility managers. From classic sturdy wood treatment tables to steel and cabinet models, Clinton tables are built for ease of use, durable performance and patient comfort.

Features:

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and raises up to 30° from flat position
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- Optional side rails & casters available
- 34" wide model available
- 550 lbs load capacity under normal use

Options:

- Power Backrest
- Headrest/Power Tables
- Safety Straps
- 2-Hook IV Pole & Mounting Bracket
- Locking Casters
- Power Footrest
- Safety Rails



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.