



Clinton Multi-Use, Ultrasound Table with Stirrups - Fully Refurbished

The Clinton Multi-Use Ultrasound Table with Stirrups is the perfect option for facilities looking for an economical way to complete basic imaging procedures in-house. Using the handheld remote control, clinicians can easily adjust the table's height and Trendelenburg / Reverse Trendelenburg. The backrest and headrest can also be adjusted for exact patient positioning. Constructed with an open base design, the Clinton Multi-Use Ultrasound Table with Stirrups allows optimal patient access. The dual lift system is for extra stability, and the all-welded base and premium mattress upholstery ensure the table will last through years of use.

Features:

- Three-piece padded top features pneumatically adjustable backrest and stirrups.
- Manual footrest can be adjusted from either side of the table and raises up to 30 degrees.
- Handheld remote control adjusts table height and Trendelenburg positioning.
- Table height ranges between 21.5 to 37.5 inches.
- Unobstructed base design
- Dual column lift system
- All-welded base
- 550-pound weight capacity under normal use.
- Comes standard with handheld remote control and paper dispenser.
- Padded headrest, arm boards, safety rails, wider mattress top and locking casters are available as additional options
- Paper dispenser included



Specs:

- Dimensions: 76"(L) x 27"(W) x 21.5" – 37.5"(H)
- Trendelenburg / Reverse Trendelenburg plus or minus 15 degrees.
- Backrest: 38"(L) x 27"(W)
- Footrest: 12"(L) x 18.5"(W)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.