



auxoMEDICAL

ConMed Sabre 180 - Fully Refurbished

This medium-powered generator delivers impressive electrosurgical performance in all modes including monopolar cut and coag and bipolar coag. Additionally, it is equipped with ARM Contact Quality, the most sensitive pad monitoring system available. The Sabre 180 is ideal for use in surgery centers and office clinics during minor, general, and laparoscopic procedures. Perfect for head & neck surgery, leep procedures, cardiac cath labs, laparoscopy carts, minor surgery centers, endoscopy, interventional radiology, and physician's offices.

Features:

- Digital Knob Power Controls
- Last Used settings on Power Up
- Aspen Return Monitor provides an extra margin of safety against possible patient burns at the return site
- The Resistance Indicator displays the resistance of the patient measured between the conductive contacts of the return electrode
- Isolated Output Circuitry minimizes the probability of alternate ground site burns
- Only one output is active at a time for maximum patient and operator safety
- Power Adjustable While Activated
- Versatile, Effective Monopolar Cutting
- Circuit Redundancy
- Microprocessor Control
- Fulguration/Coagulation
- Improved Bipolar
- Simplified Controls
- Rugged, high-efficiency circuitry
- Circuit Protection
- Mounting Flexibility
- Video Tape
- Dimensions: 6" (H) x 12" (W) x 16" (D)
- Weight: 20 lbs. (9.1 kg)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.