



## Conmed Hyfrecator 2000- Fully Refurbished

The Hyfrecator 2000 features 3 programmable power settings. It has a remote control handpiece to change power settings, maintain the sterile field, and decrease procedural time. Dual microprocessors: one to control the function and output of the Hyfrecator; the second to terminate activation if a hazard or safety issue is detected.

### Features:

- Power up/down handswitching pencil w/ 10ft. cord
- Wall Mount Kit
- Box of 50 Electrolase Disposable, non sterile blunt tips
- Box of 50 Electrolase Disposable, non sterile sharp tips
- Sterile and non-sterile sheaths
- Operator's manual, warranty card
- Two Videos: "Common Dermatologic Procedures"
- Three programmable power settings.
- Remote control hand-piece to change power settings, maintain the sterile field, and decrease procedural time.
- Dual Microprocessors: One to control the function and output of the Conmed Hyfrecator; the second to terminate activation if a hazard or safety issue is detected.



### Specs:

- Depth: 4.0" (102 mm).
- Width: 8.75" (222 mm).
- Height: 7.5" (190 mm).
- Weight: 6 lbs (2.7 kg).
- Output Power: Low 0 - 20 watts (1/10 watt increments up to 10 watts).
- High: 0 - 35 watts.
- Bipolar: 0 - 35 watts.
- Voltage: 110 V.

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.