



Covidien Valleylab Force Triad- Fully Refurbished

The ForceTriad™ energy platform is a full-featured electrosurgical system that provides electrosurgical cutting and coagulation, bipolar functionality, and vessel sealing in a single generator. The ForceTriad™ energy platform is the industry's only full-featured energy platform with remote software upgrade capabilities.

Features:

- TissueFect™ sensing technology is a Covidien control system designed to precisely manage energy delivery
- Enhanced LigaSure™ tissue fusion technology permanently fuses vessels up to and including 7 mm, lymphatics, pulmonary vasculature and tissue bundles
- An optional bipolar resection with saline software application
- Automatic instrument recognition
- All-in-one unit is compatible with conventional electrosurgical instruments and all current and new LigaSure™ instruments
- Fusion cycles are faster than the original LigaSure™ vessel sealing generator
- Valleylab™ mode provides a unique combination of monopolar hemostatis and dissection
- PER of 98 in Cut
- Connects with RapidVac™ smoke evacuator
- Adaptive REM™ system
- Offers the surgeon a wide range of options
- Easy-to-use touch screen set up
- Reduces demands on staff during procedures
- Width: 18 in. (45.7 cm)
- Depth: 20 in. (50.8 cm)
- Height: 10 in. (25.4 cm)
- Weight: 30 lb. (13.6 kg)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.