



auxoMEDICAL

Criticare nGenuity CO2 - Fully Refurbished

The Criticare nGenuity 8100E Series patient monitor is ideal for monitoring pulse oximetry (SpO₂), non-invasive blood pressure, heart rate, temperature, respiration and ECG during conscious sedation or post-procedure recovery. Dedicated function keys allow for quick setup and easy use under demanding time constraints. Full compatibility with popular Electronic Medical Records (EMR) systems makes it easy for any facility or office to transfer patient data. With a compact, lightweight design and convenient carry handle, the nGenuity 8100E monitor can be wall mounted or quickly attached to a rolling stand for easy transportation within any medical setting. For increased portability, the sealed battery can be recharged with a standard AC plug or powered directly from AC power to eliminate the need for a special adaptor.

Features:

- Monitors pulse oximetry (SpO₂), non-invasive blood pressure, heart rate, temperature, respiration and ECG.
- Easy to read 10.4 inch color display with mini-trends display of condensed tabular vital signs data.
- Large, membrane style keys
- Simplified user interface
- Compact and lightweight design
- Can be carried, wall mounted or quickly attached to a rolling stand.
- Compatible with popular Electronic Medical Records (EMR) systems
- Proprietary ComfortCuff NIBP by Criticare which measures on inflation for a faster, more comfortable reading.
- Monitors End-tidal CO₂ and inspired CO₂ accurately at high respiration rates.
- Weight: 14 lbs; 6.4 kg
- Dimensions: 11 in (H) x 13 in (W) x 10.3 (D)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.