



Curlin 4000 CMS - Fully Refurbished

The Curlin Medical 4000 CMS is a unique, technologically advanced, multi-therapy ambulatory electronic infusion pump that continuously monitors and displays "Line Pressure" to help you detect potential catheter patency problems such as occlusions, obstructions and pinched off catheters. Utilizing non-DEHP, all PVC administration sets, the 4000 CMS has the ability, at low rates, to deliver a resolution of flow that is up to 25 times better than the leading ambulatory pumps. Additionally, the 4000 CMS works synergistically with a PC and/or Palm PDA to improve the outcome, and reduce the risk of programming errors.

Features:

- Modes of Operation: Continuous, PCA with IV, Epidural, or Sub-Q Delivery, TPN with ramping, Intermittent, Variable with 24 specified doses
- Accuracy: $\pm 6\%$
- Power Sources: 2 alkaline batteries; AC and external rechargeable battery pack.
- Battery Life:
 - Nominal 10 hrs at 400 ml/hr
 - Nominal 30 hrs at 125 ml/hr
 - Nominal 85 hrs at 2 ml/hr
- Alarms: Infusion Complete, Air-in-line, Occlusion Upstream, Occlusion Downstream, Door Open, Set Not Installed, Unattended Pump, Replace Set, High Upstream pressure, Empty Battery



Specs:

- Height: 5 inches
- Width: 4 inches
- Depth: 2.5 inches
- Weight: 17.5 ounces
- Flow Rates: 0.1 ml to 400 ml/hr
- Volume Limits: 1.0 ml to 9999 ml
- KVO Rates: 0.0 to 10 ml/hr; default to 0.1 ml/hr

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.