



Datascope Passport V Patient Monitor — Fully Refurbished

Mindray's Datascope Passport V is built on a solid foundation of experience and knowledge. It follows a proven approach to acquisition and presentation of patient data, thereby delivering impressive reliability and performance. It delivers just the right combination of simplicity and sophistication, offering a comprehensive feature-set in an uncomplicated, easy-to-learn manner.

Features:

- Includes 5-lead ECG, Masimo Set® SpO2, NIBP, respiration, temperature, and a 3-trace integral recorder
- 12.1" color LCD
- 3 to 8 waveforms with automatic configuration
- Dual invasive pressures
- Lead-selectable respiration
- Continuous temperature
- Drug calculations
- 96 hours of graphic and list trends
- 128 minutes of Oxy-CRG trends
- Remote view with multi-bed alarm notification
- ECG and IBP analog output
- Device connectivity via DIAP and ethernet
- Intra-aortic balloon pump interface
- Integral 3-trace recorder



Additional Options: Arrhythmia Analysis, Nellcor® Oximax®, Microstream EtCO2, dual invasive blood pressures

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.