



Datascope Passport V with CO2 - Fully Refurbished

Passport V is built on a solid foundation of experience and knowledge. It follows a proven approach to acquisition and presentation of patient data, thereby delivering impressive reliability and performance. Enhanced capability, more optional features and extraordinary value make Passport V a compelling alternative for any healthcare environment requiring continuous physiological monitoring. Like other products in its line, Passport V utilizes a straightforward, familiar interface so that the user's experience is a positive one. It delivers just the right combination of simplicity and sophistication, offering a comprehensive feature-set in an uncomplicated, easy-to-learn manner

Features:

- 12.1" color LCD
- 3 to 8 waveforms with automatic configuration
- Enhanced numerics mode
- 3 or 5-lead ECG with or without ESU noise rejection
- Non-invasive blood pressure
- Dual invasive pressures
- Masimo SET[®] SpO₂
- Lead-selectable respiration
- Continuous temperature
- Drug calculations
- 96 hours of graphic and list trends
- 128 minutes of Oxy-CRG trends
- Remote view with multi-bed alarm notification
- ECG and IBP analog output
- Device connectivity via DIAP and ethernet
- Intra-aortic balloon pump interface
- Nurse call interface
- Defibrillator interface
- Integral carrying handle
- Front-panel alarm light
- Lithium ion battery support
- Color-coded patient connectors
- Connectivity to secondary display
- Integral 3-trace recorder



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.