



Datascope Passport XG - Fully Refurbished

The Datascope Passport XG Patient Monitor is a refurbished, multi-parameter monitor perfect for the operating or recovery setting. It is a versatile vital signs monitor specifically designed for surgery, bedside and transport applications. The parameters include ECG, Respiration, Automatic Non-Invasive Blood Pressure, Pulse Oximetry and Recorder. This unit has user-friendly menu buttons, a large screen for easy viewing, battery back-up, programmable limits and alarms and a thermal recorder. The Passport also has an easy-to-understand menu and single dedicated function keys which make it fast and easy-to-use. Includes a patient cable, power cord/AC adapter, new battery, finger sensor and standard adult blood pressure cuff.

Features:

- 3 or 5 Lead ECG
- Printer
- Non-invasive blood pressure
- Pulse oximetry
- Battery back up for transport
- SpO2
- Thermal Recorder

Specs:

- Dimensions: 12.5 x 7 x 9 inches
- Screen Size: 10.4 inches
- Batter: Type: 2 sealed lead acid batteries
- Recharge time: 16 hours
- Weight 12 Lb



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.