



Datex-Ohmeda S5 Compact - Fully Refurbished

A complete monitor for anesthesia, the Datex-Ohmeda S/5 anesthesia monitor is a multi-parameter monitor. The Datex-Ohmeda S/5 anesthesia monitor features a large screen for easy viewing in the operating room. It is an ideal monitor for any setting, including the hospital, surgery center, small OR and recovery. The Datex-Ohmeda S/5 anesthesia monitor is a top of the line unit from the leader in anesthesia monitoring.

Features:

- Small footprint bedside monitor with extra large 12.1 inch display
- Large waveforms and numerics for easy viewing from long distances and multiple viewing angles
- ECG 5-lead with ST analysis
- Impedance respiration
- Pulse oximetry (SpO2)
- NIBP
- Temperature
- Airway Gases
- CO2 (N2O), Patient O2, N2O
- 4 display options with integrated monitor controls: 15" or 17" high resolution color video display at 10.4" or 15" flat panel color LCD display.
- User configurable display layout and colors.
- Anesthesia dedicated Alarm System.
- Automatic alarm limits functions.
- 6 user configurable modes and special Cardiac Bypass Mode.
- 24 hour Numerical and Graphical trends.
- Split screen capability to show ST, Spirometry Loop, EEG, EP, or Mini Trends permanently.
- 4 user-configurable Graphical and Numerical trends.
- Waveform snapshots to store and recall waveform samples.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.