



Drager Fabius GS Anaesthesia Machine- Fully Refurbished

In an environment increasingly focused on cost issues, the anesthesia workplace Fabius GS offers an unprecedented price-performance ratio that will reduce your total cost of ownership. The ergonomic design of the Fabius GS facilitates efficient use and contributes to a positive, productive anesthesia workplace. Additionally, its modular concept allows you to configure precisely the workstation you require.

Features:

- Proven ventilation versatility
- Electronically gas-flow measurement
- Compact, convenient breathing system
- Intelligent cable management
- Universal mounting platform
- ClinicalVision

Specifications:

- Weight: 224 lbs.
- Dimensions: 35.2" 51.2" x 32.3"
- Breathing frequency: 4 to 60 1/min
- Pressure limiting (Pmax) 15 - 70 cmH₂O
- Tidal Volume
 - (Vt) 20 - 1400 mL in Volume Control
 - 20 - 1100 mL in SIMV/P
- Trigger 2 - 15 L/min
- Inspiratory pressure (P_{insp}): PEEP + 5 to 65 cmH₂O
- Pressure Support Level (Δ PPS): PEEP + 3 to 20 cmH₂O



Additional Options: Additional Drager Omega or Omega-S Monitors, Drager Vamos or SCIO Monitor

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.