



Drager Narkomed MRI-2 - Fully Refurbished

The Drager Narkomed MRI-2 is a unique anesthesia system designed with an integrated ventilator and VPO monitor that is capable of operating in shielded or unshielded MRI environments of up to 3.0 Tesla — shortening the distance between the patient and the machine. Even when placed in direct contact with the MRI, the Narkomed MRI-2 works dependably and produces no detectable image distortion. It does this while providing a complete anesthesia for a wide range of clinical operations.

Features:

- Integrated monitor measures and displays airway pressure, volume and oxygen data with airway pressure waveform.
- Customizable alarms for each unique situation.
- Standard Narkomed patient suction system with regular or reusable canister designed for MRI use.
- 3-hour battery backup
- Accurate ultrasonic flow sensor has no moving parts for easy cleaning.
- Time-cycled, pressure-limited AV2+ ventilator equipped with pneumatic circuitry and independent breathing rate and tidal volume controls.
- Vapor 2000 needs no annual recalibration, compensates for variances in pressure, temperature and flow.
- Open or closed waste scavenging system.
- Dual canister absorber
- Compatible with halothane, isoflurane, enflurane and sevoflurane anesthetic agents.

Specs:

- Dimensions: 31.5"(W) x 53"(H) x 27"(D)
- Weight: 230 lbs (104.3 kg)
- Weight of Battery Charger: 3 lbs (1.36 kg)
- Reserve Power Time: ≤ 3 hours
- Charging Time: ≤ 16 hours



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.