



ENOCHS 356 Exam Table - Fully Refurbished

ENOCHS is known as a manufacturer of a quality line of medical examination tables that proves a competitive alternative to others. The ENOCHS 356 is designed with an extra long pull out footstep that features a steel stabilizing leg designed to accommodate a 600 pound patient load capacity. The new field reversible side drawers allow the table to have the flexibility to be placed in any room configuration. The spring assisted counter balance lift mechanism allows for quick and convenient adjustment of the couch back to an infinite number of positions.

Features:

- 600 lbs. (272 kg) Patient Load
- Field Reversible Side Drawers
- Dual Electrical Outlets features locking safety door
- Instrument Warming Drawer warms drawer contents to approx. 98°
- 11° Pelvic Tilt with release levers on both sides of table
- Spring Assisted Counter Balance for back adjustment
- Extra Long, Pull-Out Locking Step with rubber tread
- Easy to clean One-Piece Upholstered Top designed without cross-seams, decorative stitching or embossing
- All welded Steel Base treated for resistance to rust
- Irrigation Basin moves independent of Leg Extension for fast transition
- Easy to adjust Stirrups with length and lateral adjustment
- Five spacious Drawers with 3.32 cubic feet of Storage
- Certified by UL to Standard 60601-1 CAN/CSA C22.3 No. 601-1
- 4-YEAR Warranty
- Back lift mechanism: horizontal to 80°
- Leg extension: 18" from end of seat section
- Paper Roll Holder: Accommodates 18" or 21" paper widths



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.