



ES Series Scrub Sink Single Bay - Fully Refurbished

The ES series scrub sinks were developed to be the answer to space and economical constraints. The ES series sinks have a smaller outside dimension but still have a deep basin to allow for a proper scrub. Made with 14 gauge 304 stainless steel that is polished and ground to a #4 finish for easy cleaning.

Features:

- ES25 has 25" wide sink
- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve (one per station, mounted under sink)
- Knee operated water (optional with IR sinks)
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Infrared sinks are equipped with a 24 volt plug in transformer
- Z-brackets (See specifications on rear panel for wall support)
- 304 stainless steel polished and ground to a #4 finish
- Net Weight of sink: 75 Lb
- Unit footprint 25"W x 20.5"D x 43"H
- Outgoing water temperature range 90-115F, factory max 115F



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.