



auxoMEDICAL

Edan DUS 60 VET Ultrasound - New

The Edan DUS 60 VET is a popular high performance, affordable portable digital ultrasound. It has now been improved with exclusive patented technology that drastically enhances the image quality for the same low price. Specifically designed for dog, cat, equine, bovine and ovine assessment, the DUS 60 VET provides robust diagnostic information. It offers a variety of multi-frequency transducers, abundant measurement and calculation software packages and Pulse Wave Doppler.

Features:

- Phase inversion harmonic imaging improves contrast resolution
- Multi-beam technology captures more detail
- Speckle reduction removes the "noise" to produce a clearer image
- Phase inversion harmonic imaging improves contrast resolution
- Synthetic aperture imaging technology enhances the signal to create better images
- Digital Beam Former
- Multi-frequency transducer series (2.0MHz - 10.0MHz)
- 12.1" high-resolution LCD monitor
- Two transducer connectors
- Two USB ports support storage, transferring and printing
- 504MB permanent image storage capacity
- 256-frame cine loop memory
- THI (Tissue Harmonic Imaging)
- TSI (Tissue Specific Imaging)
- IP (Image Processing) function
- Full QWERTY backlit keyboard
- Accessory ready: Printer or secondary monitor
- DICOM 3.0 (optional)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.