



auxoMEDICAL

Edan Elite V8 - Fully Refurbished

The Edan elite V8 Modular Patient Monitor is engineered specially for high-acuity divisions. The elite V8 dedicates to bringing high-quality healthcare to intensive cares and anesthesia monitoring, integrating world-leading technologies into one unit. The plug-and-play modular design of the elite V8 brings flexibilities to critical cares. Meantime, it offers an easy click-and-switch function to switch between modules with its intuitive MeasureSet Management interface.

Features:

- 17" high resolution color TFT-LCD screen
- Pacemaker detection, electrosurgical interference proof
- Unique iSeap algorithm specially optimized for arrhythmia patients
- Dual-mode anti-interference pulse oximetry
- VGA output and analog output
- Nurse call function & LAN connectivity
- 3/5-lead ECG, SpO₂, RESP
- NIBP, 2-TEMP, PR
- Rechargeable battery
- Optional features: 12-lead ECG, Nellcor OxiMax SpO₂, 8-IBP, C.O., Respiration CO₂ (Mainstream & Sidestream), PHASEIN Anesthesia Gas/O₂, Barcode Scanner, Touch Screen
- Defibrillation synchronization output
- Ethernet interface
- Bi-directional communications with central monitoring system
- Waveform color adjustable
- Customized alarm setting
- 720 seconds full-disclosure waveforms review.
- 150-hour graphic and tabular trends of all parameters.
- Full touch-screen enables intuitive operation by 'touching' on the specific parameters or waveform in real time
- Large font display



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.