




Edan IM80 - Fully Refurbished

The iM80 is dedicated to high-accuracy and user-friendly healthcare solutions especially for high-acuity divisions. With the latest medical and information technologies, the iM80 brings you top-quality healthcare solutions with friendly user interface. The iM80 employs various leading technologies from EDAN and other industry leaders. Its different parameter combinations can fulfill your different needs no matter it's for cardiac patients or during surgeries.

Features:

- 15" high resolution color TFT display with maximum 13 waveforms
- Touch-screen option
- Standard configuration includes the multi-purpose navigation knob
- Pacemaker detection, electrosurgical interference proof
- Defibrillation protection and defibrillation synchronization
- SD card slot enables memory extension for large data storage
- USB and serial ports enables abundant future upgrades
- VGA output and Analog output
- OxyCRG available to judge the respiration and circulation function for neonates
- Large font display
- Arrhythmia and S-T segment analysis
- Built-in rechargeable Li-ion (up to 6 hours)
- Nurse Call function and bi-directional communication with MFM-CMS central station
- Thermal Recorder
- HL7 support
- Weight: 15.4 LBS
- Dimensions: 370mm(L) x 175mm(W) x 320mm(H)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.