



Edan M3 - Fully Refurbished

The M3 has made its mark in the out-patient department, doctor's office, and general wards for its accuracy, durability, and cost-effectiveness of non-invasive SpO2, NIBP, and fast temp monitoring. Its affordable price and multi-parameter functionality address vital signs monitoring needs. The monitor is intended to be used only under regular supervision of clinical personnel. It is equipped with alarms that indicate system faults (such as loose or defective electrodes), physiologic parameters that have exceeded the limits set by the operator, or both.

Features:

- 5.7" high resolution, full color TFT Screen display for easy reading.
- Lightweight, portable design, and user-friendly interface for easy operation.
- SpO2, Pulse Rate, NIBP, and fast Temperature measurement.
- Nellcor OxiMax SpO2.
- Auto/Manual/Continuous/Average BP mode.
- Fast temp thermometer to obtain temperature in approximately 15 seconds.
- Display numeric and waveform information simultaneously.
- Real-time waveform printing.
- Backlight control and standby mode for power saving.
- Nurse call.
- Powerful storage capacity.
- Built-in lithium-ion battery for eight hours of continuous work.
- Suitable for adult, pediatric, and neonatal patients.
- PR measurement (from SpO2/NIBP).
- Trend table review and record.
- Trend graph review and record.
- USB data storage and review.
- Barcode scanner support.
- Wired and wireless network capability.
- Size: 200.8mm (L) x 241mm (H) x 189 mm (D)
- Weight 3 kg



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.