



## Edan SE-1010, PC-based - Fully Refurbished

PC-based ECG SE-1010 provides you an affordable and portable 12-channel ECG solution with powerful data analysis functions. Optional advanced analysis functions secure clinical accuracy and maximum decision-making efficiency.

### Features:

- Simultaneous 12 Leads acquisition handheld design
- Real-time ECG waveforms display on PC screen
- Auto/Manual measurement and interpretation
- WORD, PDF and JPG format reports transferred via E-mail
- Optional advanced analysis functions: Timed vectorcardiogram, VCG(3D), SAECG, QT interval dispersion, HRT, HFECG, Frequency Spectrum cardiogram
- 30s-900s HRV (heart rate variability) sample and analysis
- Stress test and exercise blood pressure monitor (option)
- Stress ECG Software
- Display Content: System name, Patient ID, Patient name, Hear rate, Display mode, Printing mode, Speed, Gain, Lowpass Filter, Hint information, ECG waves



### Specs:

- Dimensions: 5.8" x 3.9" x 1.6"
- Leads Mode: 12 standard leads/cabrera leads
- Acquisition Mode: Simultaneous 12 leads acquisition
- Sampling Rate: 1,000 sec/channel
- A/D Resolution: 24 bits
- Time Constant:  $\geq 3.2$  s
- Frequency Response: 0.05Hz - 150Hz

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.