



# auxoMEDICAL

## Edan SE-1200 ECG- Fully Refurbished

Excellent cost performance with built-in high-resolution thermal printer. You can record 12 ECG traces on A4 size and US Letter-size paper. 5.7" color touch-screen, Alphanumeric keyboard and one touch operation, and Sleep mode to save energy and extend LCD life.

Automatic measurement and interpretation tested with authoritative CSE database along with Real-time waveform freezing. 120 seconds ECG waveform reviewable, Two-step exercise test with periodic recording, and Built-in rechargeable Li-ion battery.

### Features:

- 5.7 inch screen for optimal viewing
- Alphanumeric keyboard and one touch operation
- Sleep mode to save energy and extend LCD life
- Automatic measurement and interpretation tested with authoritative CSE database
- Real-time waveform freezing
- 120 seconds ECG waveform reviewable
- Two-step exercise test with periodic recording
- Built-in rechargeable Li-ion battery
- 100 ECGs in internal memory
- Optional: PC-based data management software (Smart ECG Viewer)
- One-touch operation
- Enhanced weak signal detection
- Comprehensive filters & anti-noise technology
- Advanced SEMIP algorithm for automatic measurements and interpretation
- Built-in high resolution thermal printer with external USB printer (optional) support
- Extended printing
- Large internal memory which can be enlarged by USB flash disk
- Built-in data management with reports review
- LAN/Wifi(optional) connections
- Dimensions: 420mm×330mm×105mm (16.5in×13.0in×4.1in)
- Weight: Approx. 5.0kg (11.0 lbs)



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.