



auxoMEDICAL

Edan SE-601 Series ECG - Fully Refurbished

The Vista Basic was designed to ensure simplicity in all aspects of its use: from set installation, to programming, to special mode use and alarm correction. Based on a design that has sold over 100,000 units worldwide, the Vista Basic provides high reliability packaged in one of the easiest pumps to use. As a workhorse in the outpatient market, it offers a rugged and reliable design, all relevant therapies such as continuous, piggyback and program mode, and a selection of cost-effective, straight line gravity sets.

Features:

- Compact, lightweight and portable design
- Vivid waveform display on 5.7" LCD screen
- Alphanumeric keyboard and one touch operation
- Pacemaker detection for the requirements of ANSI/AAMI ECII
- Automatic measurements and interpretation tested by authoritative CSE/MIT/AHA database
- Heart Rate Variability (HRV) analysis
- Internal thermal printer and external printer via USB port
- PC-based data management software (Smart ECG Viewer)
- Data transmitting to PC via LAN/WIFI
- One-touch operation
- 12 leads simultaneous acquisition
- Enhanced weak signal detection
- Comprehensive filters & anti-noise technology Models
- Dimensions: 288mm x 210mm x 70mm
- Input Circuit: Defibrillation protection built-in
- Time Constant: >3.2s
- Record Paper: 110mm x 140mm* 144P (Z-fold)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.