



auxoMEDICAL

Erbe ICC 200 - Fully Refurbished

The ERBIE ICC 200, a comfortable unit for all those who need no more than 200 watts of output, with CUT CONTROL for sharply defined, reproducible cutting quality. Optional ENDO CUT gives the best in cut control and fractionation in polypectomy or papillotomy. A comprehensive range of coagulation functions

Features:

- ENDO CUT™ for use with snare wires and sphincterotomes
- Power on Demand – As much as necessary, as little as possible
- Forced Coag™
- Soft Coag™
- Expandable for Argon use
- Compatible with:
 - Snare Wires
 - Sphincterotomes
 - Hemostasis Probes
 - Hot Biopsy Forceps



Specs:

- Dimensions (HxWxD) cm, (In)
(Display) 28 x 15.2 x 36.8 (11 x 5.9 x 14.5)
- Weight: 8kg (17.6 lb)
- Generator: Solid-state
- Frequency, khz: 350
- Power Neede, VAC (Energy): 80-140 (110/230)

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.