



## FHC 7100 Mobilecare - Fully Refurbished

The FHC7100 Mobilecare Transport Hospital Stretcher is designed to provide deluxe transportation features without the deluxe price tag. With 8" Tente total/directional locking casters and brake steer pedal at every corner, this transporting stretcher provides easy transportation of all patients up to 550 lbs. The 7100 stretcher comes with I.V. pole sockets at every corner and fold forward side rails with built in push handles for convenient transporting needs.

A O2 storage underneath the hospital stretcher make it simple to keep patient oxygen stored while moving. With our FHC7100 Mobilecare Transport Hospital Stretcher, moving sick or injured patients around couldn't be any easier or simpler.

### Features:

- Pneumatic back adjustment
- Height adjustment pedals on both sides from 23" -36" high
- Trendelenburg & Reverse Trendelenburg adjustment pedals on both sides
- Fold down siderails
- Built-in O2 cylinder storage
- Central locking casters
- Brake/Steer pedal at each corner
- IV pole
- New stretcher pad
- Patient Surface: 75" L x 26" W
- Patient Weight Capacity: 550 lbs.



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.