



## GE Aespire 7100 Anesthesia Machine - Fully Refurbished

The Aespire 7100 is a compact, integrated workstation that helps simplify the daily demands of your work environment. With volume ventilation mode as standard, the Aespire 7100 can be enhanced by several software options to help you optimize clinical performance in an economic package.

### Features:

- Enhanced monitor integration capabilities with out Cardiocap 5 and Dash series of monitors
- Large color ventilator display with color waveform and alarm message indicator
- Lightweight and compact
- Optional integrated auxiliary O2 flowmeter and suction control
- Advanced Breathing System (ABS)
- Volume and pressure modes with electronic PEEP
- Improved low flow/reduced life cycle costs
- Equipped with Volume Control ventilation, integrated volume and oxygen monitoring, electronic PEEP and alarm management

### Specifications:

- Dimensions: 28.3" x 28.7" x 52.9"
- Weight: 238 lbs.
- Tidal volume range: 45 to 1500 mL
- Pressure range: 5 to 50 cm H2O
- Rate: 4 to 65 breaths per minute

**Additional Options:** Pressure Mode, Pressure Waveform, Tidal Volume Compensation



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.