



## GE Aespire with 7900 Smart Vent - Fully Refurbished

The Aestiva/5 7900 anesthesia system combines the premiere features of the Aestiva/5 anesthesia system with GE's newest anesthesia ventilation platform - the 7900 Ventilator. The 7900 Ventilator has been specifically engineered to allow for ventilation of a wide range of patients making it an ideal solution for clinicians that are looking for a high performance ventilator at a cost-effective price. The design of the 7900 ventilator is built around the proven SmartVent technology, incorporating similar user interface, menus and screen views. Its design helps to deliver the level of performance and patient safety which customers have come to expect from our products.

### Features:

- 7900 Vent w/4.8 SW. Press and Volume control Vent Modes
- PSV Pro and SIMV Vent Modes Included (4 modes total)
- 3 gas delivery system
- Optional suction control
- Lightweight and compact, easy to maneuver
- One step bag/vent switch turns the ventilator on/off
- Smart alarms direct user to specific problems and affected parameters
- Improved low flow/reduced life cycle costs
- Dual air flow tubes standard for higher resolution of low flows
- EZ Change Soda Sorb Canister system
- Enhanced hose and cable management helps reduce disconnects and clutter.
- Large work surface offers more room for supplies and chart.
- Large drawers provide exceptional storage.
- Removable drawer insert conveniently stores small items



### Specs:

- Dimensions: 53.5"(H) x 28.5 in"(W) x 28.9 in (D)
- Weight: 136 kg (300 lb)
- Top shelf weight limit: 35 kg (75 lb)
- Minimum O2 flow of 50 mL

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.