



GE Datex-Ohmeda S/5 Aespire 7900 - Fully Refurbished

The S/5 Aespire® was developed using GE's unique approach to anesthesia workstation development – close and continuous collaboration with clinicians. It combines the proven performance of the 7100 Ventilator or 7900 Ventilator with Datex-Ohmeda® monitoring technology into an affordable system solution designed for clinicians throughout the world. The S/5 Aespire® utilizes the Advanced Breathing System, the ABS™. Its design, based on your needs and expertise.

Features:

- Optional suction control eliminates need for stand alone suction
- Optional auxiliary O2 supplies a pure oxygen source for use with a nasal cannula
- LED light strip provides bi-level work surface illumination
- Auxiliary Common Gas Outlet facilitates low pressure leak test
- Precision, reliable pneumatic gas flow mixing
- Volume Control
- Pressure Control
- Synchronized Intermittent Mandatory Ventilation (SIMV) with Pressure Support
- PSVPro® (Pressure Support with Apnea backup mode)
- Power up is quick and easy
- Individual flow controls
- Enhanced hose and cable management
- Larger work surface
- Two large drawers
- Easy to clean, fully autoclavable, latex-free
- Integrated electronics
- Faster response - ideal for low flow anesthesia
- Multi Absorber canister facilitates fast, easy removal and replacement
- One step mechanical bag to vent switch turns ventilator on/off requires no electronics



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.