



GE Healthcare Aestiva/5 - Fully Refurbished

A used or refurbished GE Aestiva®5 is reliable, fully featured, and can be a very cost effective way to deliver anesthesia. The GE Datex Ohmeda Aestiva®/5 with SmartVent offers you some of the capabilities of a respiratory ventilator, but saves the cost of bringing a separate ventilator into the operating room. Its flexible design and straight forward operations enhance the SmartVent's capabilities to give a higher level of control. The Aestiva®/5 configuration gives you the ability to use your current monitors and data management systems, helping you stay within your budget while achieving physical integration in one unit.

Features:

- Volume Mode, Pressure Control Mode, Pressure Support (PSVPro®), Synchronized Intermittent Mandatory Ventilation (SIMV), electronic PEEP
- Tidal volume compensation
- One motion from mechanical to manual mode
- Two key presses to total standby: end case
- Cardiac bypass case mode
- Lower overall height
- User configurable drawers/shelving
- Eight machine hoses/cables integrated
- "No tools" disassembly of components
- Autoclavable and latex-free
- Responsive location of common gas outlet
- Fresh gas flow compensation—automatically
- Smooth, faster acting fresh gas flow control
- Minimum O2 flow of 50 mL
- Dual air flow tube for low flow
- Two scheduled maintenance checks per year
- 2 Vaporizer configuration Dimensions: 53.4" x 29.5"(W) x 32.7"(D)
- 2 Vaporizer configuration weight: 136 kg/300 lb
- 3 Vaporizer configuration Dimensions: 53.4"(H) x 36.6"(W) x 32.7"(D)
- 3 Vaporizer configuration weight: 154 kg/340 lb



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.