



# auxoMEDICAL

## GE Logiq e - Fully Refurbished

The GE Logiq e is one of the most popular portable ultrasound systems for sale. It holds this honor for a few reasons: Price, Performance, and Capability. Also benefitting the GE Logiq e is that it's extensively customizable. It allows you to adjust all aspects of the image quality, annotations, presets, programmable buttons, etc. For these reasons it remains a very popular portable ultrasound solution.

### Features:

- 15" LCD
- Versatile and shared-services with CW Doppler
- CrossXBeam compound imaging (as option)
- SRI Speckle Reduction imaging (as option)
- Excellent for MSK, orthopedic, vascular, venous, and needle guided injections
- Includes full set of applications and calculations
- Coded Harmonic Imaging
- More customizable than many of its competitors
- Trackball allows for easier measurements than trackpad systems
- 2D, M-Mode
- PW/CW/Color/Power Doppler

### Specifications:

- Height (console only): 61mm (2.49 in),
- 76.5mm (3.12 in) with handle.
- Weighs approximately 4.6kg (10.1 lb)
- Width: 340mm (13.88 in)
- Depth: 287mm (11/71 in) for console only, 327 mm (13.35 in) w/handle
- Voltage: 100-240 V AC
- Frequency: 50/60Hz



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.