



# auxoMEDICAL

## Genoray ZEN-7000 - Fully Refurbished

The Genoray Zen-7000 C-Arm is the solution to your full size C-Arm needs. For the pain management clinic, surgery center of hospital, this system is the superior product to add to your operating room. The Genoray Zen-7000 supports Windows 7 software for an easy-to-use interface. A 1K x 1K CCD Camera provides a clear picture every time. This C-Arm is loaded with features including dual medical grade monitors, Toshiba image intensifier, Boost, Rotating Anode, and is RIS MWL Compatible.

### Features:

- Rotating Anode
- Tri Mode 9"/6"/4.5" TOSHIBA Image Intensifier (I.I.)
- Auto Brightness Control
- Large Focus (Boost)
- 1K x 1K CCD Camera (proprietary Genoray)
- 110v requirement
- Touch Screen Operation Panel
- RIS MWL Compatible
- Fluoro Boost/Pulsed Fluoro
- Windows 7 Based
- 100,000 Image Storage
- LIH - Last Image Hold
- Dose Reduction Pulsed Fluoro
- Digital Record & Annotation
- CINE Recording & Contrast Adjustment
- USB Ports & Sharp Contrast
- CDR W, DICOM, Zoom, Negate & Edge Enhancement
- CINE 30 Frames per Second Digital Image Processing
- Minimum Radiation Output



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.

