



Hausted 578 Eye Stretcher - Fully Refurbished

The Gemini 578 is a multi-position stretcher intended for treatment, transport, or recovery of a patient. With the uniquely designed articulating head section, this unit can be used for specialty procedures which require precision head positioning, such as ophthalmology, plastic surgery, or ENT procedures.

Features:

- Low profile height adjustable articulating head piece with dual operation control knobs.
- Stainless steel 4-leg, 3-position Retracto rails with top caps that eliminate patient transfer gap.
- Dual side foot pedal hydraulic controls for height adjustment.
- Brake and steering activation from all 4 corners.
- 8" easy roll casters.
- Stainless steel patient surface.
- Seamless thermo-formed base cover with storage and built-in oxygen tank holder.
- Perimeter non-marring bumper system.
- Arm rests around Patient Headrest not included (pictured)



Specs:

- Trendelenburgs: 0-16°
- Height Range: High: 32" (81 cm). ; Low: 24" (61 cm).
- Overall Width: 29.8" (75.8 cm).
- Patient Surface Width: 24.5" (62.2 cm).
- Overall Length: 88.5" (225 cm).
- Patient Surface Length: 81.5" (207 cm).
- Backrest Degree of Movement: 0-90° (crank operated).
- Knee Flex Degree of Movement: 0-75° (crank operated).
- Retracto Rail Height: 14" (35.5 cm).
- Retracto Rail Length: 49" (125 cm).
- Max Weight Capacity (Evenly Distributed): 500 lbs (227 kg).

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.