



Hausted Horizon Hydraulic Series - Fully Refurbished

The Horizon Hydraulic Stretchers are designed to provide versatility, lasting and consistent performance, and ease of use. These easy-to-clean units can support a patient during induction, transport, treatment, and recovery functions, avoiding unnecessary patient transfers to other equipment. The Horizon Hydraulic Stretcher's versatility saves time and effort. It protects facility walls and the stretcher with our heavy-duty, shock-absorbing, non-marking bumpers. Easy-to-operate casters and brakes to help reduce safety risks to healthcare personnel.

Features:

- Zone-of-Care™ allows all functions to be performed at the patient's side
- Built-in central steering and braking systems
- Pedals located at all four corners.
- Hydraulic foot pedals on both sides
- Hands-free height adjustment and Trendelenburg positioning
- Dual-purpose accessory wells at all four corners
- Open base design holds oxygen tanks and patient belongings



Specs:

- Overall Width – 31.5" (80 cm)
- Side Rails – 15 x 75" (38 x 191 cm)
- Knee Gatch – No
- Rail Type – Airglide
- Patient Surface Dimensions – 28 x 76" (71 x 193 cm)
- 462-APA-ST

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.