



## Hausted VIC Video Imaging Chair - Fully Refurbished

Chair designed specifically for modified barium swallow studies for AP and Lateral imaging as a procedural chair for diagnosis and treatment. This unit adjusts easily from sitting to supine and provides a rotating seat section with a stable work surface during patient treatment. Transport, treatment, procedure, and recovery can be performed on the same unit to improve case turn around and stream-line patient handling.

### Features:

- Specially-Designed Back Section Accommodates the Space Confinements of Large Imaging Equipment
- Radiolucent Back Section Provides a Clear, Unobstructed Site For Imaging Used During Transport, Easy To Remove For Radiological Procedures
- Easy-To-Operate Rotating Seat Section Saves Time in Repositioning the Patient For Anterior-Posterior And Lateral Imaging Views
- Dual Side Hydraulic Foot Pedal Controls For Height Adjustment: 23.5" to 31.5"
- Pneumatic Assist Backrest: 0° to 90°
- Trendelenburg: – 10°
- Four Wheel Brake and Steer System Operable From Both Sides of Unit Feature 5 In Easy Rolling Casters
- Easy To Clean Stainless Steel Base Cover
- Polymer Coated Steel Litter Top
- Fold-Down Foot Rest (Height Adjustable)
- Accessory Wells (2 on Each Side)
- Built-In Chart Holder and Push Bar
- Swing Down – Tuck-Under Rails with Padded Armrests Maintain Patient Security
- Patient Surface: 22" x 70.5"
- Weight Capacity: 325 lbs



### Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.