



HeartSine Samaritan PAD 350P - Fully Refurbished

The HeartSine samaritan PAD is an automated external defibrillator (AED) that analyzes the heart rhythm and delivers an electrical shock to victims of Sudden Cardiac Arrest in order to restore the heart to normal rhythm. This user-friendly AED provides easy to follow visual and audio prompts to guide the rescuer from pad application to shock delivery, if required.

Features:

- Portable and lightweight (2.4 lbs)
- Resists shock and vibration and carries an IP56 Rating, the industry's highest rating against dust and moisture ingress
- Utilizes proprietary electrode technology, advanced and stable firmware, and proprietary SCOPE Biphasic Technology
- Optimized biphasic escalating waveform compensates energy, slope and envelope for patient impedance
- Easy-to understand visual and voice prompts that guide a user through the process, including CPR coaching
- Automatically runs self check each week
- Two button operation
- Evaluates patient's ECG, signal quality, electrode contact integrity and patient impedance to determine if defibrillation is required



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.