



Hill-Rom 1105 Advance Bed — Fully Refurbished

The 1105 Hill-Rom Advance bed is built to last while maintaining patient comfort. Along with a retractable frame and easy to clean plastic head and footboards. The 1105 Bed is the retractable type that moves the patient closer to the head end of the bed automatically when the head section is elevated. The bed is equipped with a night light.

Features:

- 500 lb capacity
- New Patient Mattress
- Foot-end controls
- High-low travel time: 21 seconds
- Sleep surface: 36" (w) x 82" (l)
- Overall dimensions (w/out bumpers): 36.5" (w) x 91" (l) • Bed frame clearance with 5" casters: 7.5"
- Low position: 17"
- High position: 31"
- Central pedals for braking and steering
- Patient night light
- Instant CPR release
- One year parts warranty



Options:

- Prestige walnut end grain
- Natural oak end grain

Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.