



Hill-Rom General Purpose Stretcher - Fully Refurbished

The GPS General Procedure Stretcher is ideal for transport, recovery, and minor procedures. It features hydraulic adjustable height, with controls located on both sides of the stretcher. The side rails collapse and can be tucked away when not in use, allowing zero transfer gap. All 4 wheels brake from 1 pedal located at all 4 corners of the stretcher.

Features:

- Hydraulic Hi-Low
- Pneumatic Fowler
- One Step Four Wheel
- Locking Casters
- Retractable Side Rails (with fold under feature)
- Thermaform Base
- Head Up/Down & Knee Up/Down
- Cardiac Chair
- Trendelenberg & Reverse Trendelenberg
- All wheels brake from any corner



Specs:

- 500 lb. patient capacity
- 33" x 80" Sleep Surface
- 28.5" or 33.5" Width
- Adjustable Height of 22.5" to 38.5"

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.