



Hill-Rom P8005 Transport Stretcher - Fully Refurbished

The Next Generation Steering Plus System on Hill-Rom's Transport Stretcher is easy to steer and stop while allowing for crisp cornering and control. Increased stability is provided through the dual-locking four wheel braking system, which is standard on all Hill-Rom Stretchers. Built to last, the transport stretcher's heavy-gauge, automotive-grade steel reinforced frame provides both long-lasting quality and a 700 lb. weight capacity to handle patients of varying size. The stretcher integrates well with other emergency room equipment, making it an ideal ER and transport stretcher.

Features:

- 700 lb weight capacity
- Low deck height capability
- Push handles
- Redundant-locking, four-wheel brakes
- Power washable
- Next generation Steering Plus system
- Permanent IV pole- 2 stage
- Overall length: 83" (210.8 cm)
- Overall width: 32" (81.3 cm) or 36" (91.4 cm)
- Siderail length: 47" (119.4 cm)
- Siderail height above deck: 14.5" (36.8 cm)
- Mattress size: 26" (66 cm) x 75" (190.5 cm) or 30" (76.2 cm) x 75" (190.5 cm)
- Maximum head elevation: 90°
- Maximum Trendelenburg/Reverse Trendelenburg: 18°



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.