



Hill-Rom P8050 OB-GYN Stretcher - Fully Refurbished

With its fast and easy patient positioning, Hill-Rom's OB-GYN stretcher (P8050) is ideal for use in both maternity and emergency.

Features:

- 700 lb weight capacity
- Sliding patient litter
- Next generation Steering Plus
- Dual-locking four-wheel brakes
- Zero Transfer Gap Tucking Siderails
- Power washable
- Sliding litter
- Integrated foot supports
- No-lift stowaway foot section
- Active hand brake
- Integrated IV pole transport device
- Brake/Steer pedals – head and foot
- Push handles
- Permanent IV pole – 3 stage (head only)
- 30-inch width
- Radiolucent deck
- 3-inch high density foam mattress with fire barrier
- Integrated telescoping calf supports
- Maximum weight limit: 700 lbs (317.5 kg)
- Overall length: 80.5" (204.5 cm)
- Overall width (siderails up): 35.5" (90.2 cm)
- Overall width (siderails stored): 31" (78.7 cm)
- Siderail length: 47" (119.4 cm)
- Siderail height from top of Deck: 13" (33.0 cm)
- Mattress size: 30" (76.2 cm) x 75" (190.5 cm)



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.