

Hill-Rom TotalCare Bed — Fully Refurbished

Developed in conjunction with leading multi-disciplined healthcare professionals, this revolutionary frame and surface system creates value by helping to reduce back-related injuries associated with patient ambulation, repositioning, and transfers. The TotalCare® Bed System is truly a single-product solution that maximizes human and capital assets.

Features:

- Line-of-Site Angle Indicators
- FullChair One-Button Chair Positioning Mechanism
- 12 Amp Accessory Outlet
- Graphical Caregiver Interface (GCI) Control
- Manual Foot Pedal Hydraulics
- HandsFree Emergency CPR & Trendelenburg Release Mechanism
- Drainage Bad Holders
- Point-of-Care Brake and Steer System
- Safety & Information Indicator Panel
- WallGuard Bumper System
- Permanent IV Poles (Both Sides)
- Transport Handles
- Removable Headboard for CPR Support
- OneStep Sidereal Release Mechanism
- Speaker
- Patient Controls



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.