



# auxoMEDICAL

## Hospira Lifecare PCA3 - Fully Refurbished

In conjunction with Hospira MedNet software, this next-generation LifeCare PCA system offers an unmatched combination of features for patient-controlled analgesia therapy, including, Hospira MedNet Software with customizable drug libraries, integral bar code reader for medication identification of prefilled and custom-filled vials with pharmacy-generated bar codes, multiple medication rule sets for as many as 18 clinical care area, bar-coded, prefilled medication vials, stored protocols for quick, accurate programming of standing orders and easy-to-use design based on human factors research.

### Features:

- Wireless connectivity capability that facilitates real-time infusion monitoring, including support of real-time 5-rights verification.
- Confirmation screens and stored protocols
- Intuitive, easy-to-understand programming screens.
- A built-in bar code reader to help confirm the right medication and the right concentration.
- Pre-filled, bar-coded medication vials with colored labels that help reduce errors and assure a correct match to the hospital's medication rule set.
- A clear, easy-to-read numeric keypad that allows direct input of numbers.
- Hospira Mednet Safety Software

### Specs:

Dimensions: 8"W 13"H 6"D (excluding pole clamp protrusion and power cord storage)

Weight: 10 lbs. with battery.

Battery: One sealed, lead-acid, rechargeable 8 V battery, internal to device.

Battery Life: 4 hours of continuous operation



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.