



Mac PT-1000 - Fully Refurbished

The MAC PT-1000 general transport stretcher allows for easy and comfortable patient transfer with a weight capacity of up to 550 lbs. Economical and versatile, this stretcher is an ideal solution for any department or facility where cost efficient quality is required.

Features:

- Brake steer pedals all 4 corners
- Dual pedestal lift mechanism
- Trendelenburg/reverse Trendelenburg
- Lift/Trendelenburg pedals on both sides
- Fold forward side rails with built in push handles
- 2 hook removable telescopic I.V. pole
- I.V. pole sockets at all corners
- O2 storage underneath
- Dimensions: Overall: 83" long x 31" wide
- Patient surface: 75" long x 26" wide
- Overall: 83" long x 36" wide
- Patient surface: 75" long x 31" wide
- Equipped with dual fowler assist
- Height Adjustment Range: Minimum 23" Maximum 35"
- Litter Positioning: Trend/reverse trend: +/- 18°
- Fowler position: 90°
- Weight capacity: 550 pounds
- Overall: 83" long x 31" wide



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.