



auxoMEDICAL

Medivator DSD-201LT - Fully Refurbished

Automated endoscope reprocessor features two fully-independent and asynchronous reprocessing basins for high-volume, quick endoscope turnaround. The MEDIVATORS® DSD-201LT has user-programmable reprocessing cycles, and an enhanced sensor package for continuous monitoring of all cycle parameters. System design self-disinfects internal fluid lines during every cycle, and a log printout documents all test results.

Features:

- Fast cycle time for quick endoscope turnaround
- Two basins for simultaneous or independent cycles
- Cycle parameters continuously monitored
- Leak test ensures endoscope protection throughout cycle
- Automated alcohol flush and air dry
- Audible and visual alarms for reprocessing cycle deficiency
- Assures patient and operator safety
- Cycle customization meets specific reprocessing needs
- Bacteria retentive three-stage water filtration system
- Compatible with all current models of endoscopes, including those with elevator wire channels
- Printed documentation of all reprocessing cycle data
- Self-disinfects internal fluid lines during every cycle
- For use with most FDA-cleared reusable high-level disinfectants
- Optional computerized endoscope reprocessing records, vapor management system and neutralization system

Specs:

- Dimensions: 46 H x 36 W x 21 D (inches)
- Height with lid open: 64 inches (162.5 cm)
- Weight: 400 lbs (182 kg)
- Reservoir Capacities: Disinfectant: 4 gallons (15.1 liters), Alcohol and detergent: 27oz each (800 ml)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.