



## Medstone Elite Arcoma BT Mobile Imaging Table - Fully Refurbished

Arcoma designed and manufactured the Medstone™ "Elite and Precision" mobile imaging tables for pain management, vascular, urology and lithotripter with emphasis on high quality and reliability for use with the high end mobile digital c-arm with complete 360® rotation of unencumbered 3D imaging. The Arcoma "Medstone™" series of imaging table have become the benchmark in the industry and are renowned for their reliability since their original introduction in 1999.

### Features:

- 76"x20" carbon fiber tabletop with 28" table extension
- Allows for 50" of unobstructed imaging
- 500lb weight capacity
- Hand and foot control w/ 5 table motions:
  1. Motor Driven Extended-retract travel: 16" (41cm)
  2. Motor Driven elevating travel: 31" min. 43" max. (69cm to 109cm.)
  3. Motor Driven Lateral movement: 8" (20cm)
  4. Motor Driven Trendelenberg tilt: +/- 15°
  5. Motor Driven lateral tilt: +/- 15°



### Accessories Included:

- Foot Control
- Hand Control
- Side Rails

## Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.