



Midmark 204 Exam Table - Fully Refurbished

The Ritter 204 Manual Exam Table features the latest in design, comfort and functionality in a manual exam table. It provides an enlarged footstep for easier access to make patients feel secure, three side drawers for optimal supply management, modern styling to coordinate with your office, an ergonomic top design to help patients feel comfortable, and additional weight capacity to handle heavy patients.

Features:

- Pass-thru drawer design
- The ergonomic design of the upholstered top means this top is as soft and comfortable as it looks
- The seamless upholstery eliminates areas where dirt and bacteria can accumulate
- Removable exam table top system
- Enlarged footstep allows patients to feel more secure and get on the table more confidently
- Spacious side drawers in addition to front drawers
- Increased weight capacity
- Height with upholstery: Ritter Soft Touch Top (32.13")
Ritter Seamless Top (33.13")
- Duplex, hospital grade electrical receptacle on left side of table
- The Ritter 204 Manual Exam Tables meet UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 Standards
- Multi-Density Foam Top Design
- Powder-coated paint finish
- 18-gauge steel shell
- Seamless polystyrene drawers
- Ball-bearing drawer glides
- 500 lb. weight capacity



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.